

BALTIC

manufacturer for OR and ICU with various product lines of surgical light, operating table, medical supply unit and digital OR system, providing safe, effective and cost-saving solutions for modern hospitals.

Specifications :

- 1.The size of table : 1900*650mm
- 2.Height adjustment : 710mm-930mm
- 3.Table board trend : Reversed trendelenburg \geq 14°
Trendelenburg \geq 17°
- 4.Back board adjustment : 0°- +76°
- 5.Power supply : AC220V, 50Hz
- 6.Power consumption \leq 200VA

Accessories :

1 Arm rest	1406.3000	-----	2
2 Fasten socket(square)	1402.2000	-----	2
3 Handhold controller	1418.0001	-----	1
4 Fasten socket (adjustable)	1401.1000	-----	2
5 Foot switch	1501.4000	-----	1
6 Fluid basin	1404.0000	-----	1
7 Foot pedal	1501.4000	-----	2
8 Leg support	1505.0000	-----	2
9 Hand grip	1406.4000	-----	2

SMD2 9100 VERA



BALTIC

Electric Delivery /OB/GYN Table



PRODUCT FEATURES



Other colors 

- New ergonomic design with multi-purpose: can be used as Delivery Bed, OB/GYN Operation Table or Examination Couch.
- Attractive appearance and electric motor controlled. Height, back board, sitting board position is electrically adjusted by foot switch or handhold controller.
- Assistant table (leg board) can be folded and hidden beneath the sitting board. Compact and space saving.
- Unique concave shape of base ensures an easy access to the patient during operation.
- Plug-in basin can be easily moved out for a thoroughly cleaning and disinfection.
- Robust steel structure table frame and casting iron table base with stainless steel cover. Providing a stable and safe support.
- Table surface and accessories are made of high class stainless steel. Durable and anti-rust, easy cleaning and disinfection.
- Installing and dismantling the sitting board easily and rapidly.
- Safe, reliable and low-noise.
- Seamless PU mattress, ensure excellent comfort for long time sitting or lying. Antibacterial, antiallergic, fluid

